

# LIMK2 antibody



Catalog Number: 112229

## Product name

LIMK2 antibody

## Immunogen

[Human LIMK2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LIMK2 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 70-80 kDa

## Preparation

This antibody was obtained by immunization of LIMK2 recombinant protein (Accession Number: NM\_001031801). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112229(LIMK2 Antibody) at dilution of 1:600