

# LINGO1 antibody



Catalog Number: 112240

---

## Product name

LINGO1 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

LINGO1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, rat brain tissue. Observed molecular weight by Western-blot: 83 kDa

## Preparation

This antibody was obtained by immunization of LINGO1 recombinant protein (Accession Number: NM\_032808). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112240(LINGO1 antibody) at dilution of 1:1500