

# Kv1.3 antibody

Catalog Number: 112256

---

## Product name

Kv1.3 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Kv1.3 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 52 kDa, 75 kDa

## Preparation

This antibody was obtained by immunization of Kv1.3 recombinant protein (Accession Number: NM\_002232). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112256(KCNA3 antibody) at dilution of 1:1000