

Kv1.4-Specific antibody



Catalog Number: 112257

Product name

Kv1.4-Specific antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Kv1.4-Specific Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Observed molecular weight by Western-blot: 73kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002233). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112257(KCNA4-Specific antibody) at dilution of 1:500