

LAG-3 antibody

Catalog Number: 112268

Product name

LAG-3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LAG-3 Rabbit Polyclonal antibody. Positive IP detected in mouse liver tissue. Positive WB detected in mouse liver tissue, A2780 cells, Jurkat cells. Observed molecular weight by Western-blot: 54kd, 60-65kd

Preparation

This antibody was obtained by immunization of LAG-3 recombinant protein (Accession Number: BC052589). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112268(LAG3 antibody) at dilution of 1:500



IP Result of anti-LAG3 (IP:Catalog No:112268, 4ug; Detection:Catalog No:112268 1:500) with mouse liver tissue lysate 4000ug.