

Product name

LRMP antibody

Specificity

Human; other species not tested.

Antibody description

LRMP Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human spleen tissue.

Preparation

This antibody was obtained by immunization of LRMP recombinant protein (Accession Number: NM_001366545). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IHC


Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations

Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112320(LRMP Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112320(LRMP Antibody) at dilution of 1:50 (under 40x lens)