

LRRC38 antibody



Catalog Number: 112332

Product name

LRRC38 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

LRRC38 Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM_001713942). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112332(LRRC38 Antibody) at dilution of 1:300