

# LRWD1 antibody



Catalog Number: 112342

---

**Product name**

LRWD1 antibody

**Specificity**

Human; other species not tested.

**Antibody description**

LRWD1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 85 kDa

**Preparation**

This antibody was obtained by immunization of Peptide (Accession Number: NM\_152892).  
Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG

**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Validations**

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112342(LRWD1 Antibody) at dilution of 1:600