

# LSM10 antibody



Catalog Number: 112346

## Product name

Polyclonal

LSM10 antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

LSM10 Rabbit Polyclonal antibody. Positive WB detected in human spleen tissue. Observed molecular weight by Western-blot: 14kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of LSM10 recombinant protein (Accession Number: NM\_032881). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:112346(LSM10 antibody) at dilution of 1:200