

LSM11 antibody



Catalog Number: 112347

Product name

LSM11 antibody

Specificity

Human; other species not tested.

Antibody description

LSM11 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 45-50 kDa

Preparation

This antibody was obtained by immunization of LSM11 recombinant protein (Accession Number: NM_173491). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112347(LSM11 Antibody) at dilution of 1:1000