

LSR antibody

Catalog Number: 112354

Product name

LSR antibody

Specificity

Human,Rat, Mouse; other species not tested.

Antibody description

LSR Rabbit Polyclonal antibody. Positive FC detected in HepG2 cells. Positive WB detected in rat liver tissue, L02 cells, mouse liver tissue. Observed molecular weight by Western-blot: 56 kDa, 68 kDa

Preparation

This antibody was obtained by immunization of LSR recombinant protein (Accession Number: BC004381). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



rat liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112354(LSR antibody) at dilution of 1:1000



1X10⁶ HepG2 cells were stained with .2ug LSR antibody (Catalog No:112354, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.