

LTA4H antibody



Catalog Number: 112357

Product name

LTA4H antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

LTA4H Rabbit Polyclonal antibody. Positive IHC detected in human skin tissue, human brain tissue, human heart tissue, human kidney tissue, human lung tissue, human ovary tissue, human placenta tissue, human spleen tissue, human testis tissue. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 57-70 kDa

Validations




Preparation

This antibody was obtained by immunization of LTA4H recombinant protein (Accession Number: NM_000895). Purification method: Antigen affinity purified.

Immunohistochemical of paraffin-embedded human skin using Catalog No:112357(LTA4H antibody) at dilution of 1:100 (under 10x lens)

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Immunohistochemical of paraffin-embedded human skin using Catalog No:112357(LTA4H antibody) at dilution of 1:100 (under 40x lens)

Clonality

Polyclonal

Ig Type

Rabbit IgG



Applications

ELISA, IHC, WB, IP

Dilutions

LTA4H antibody

Catalog Number: 112357



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112357(LTA4H Antibody) at dilution of 1:1000

IP Result of anti-LTA4H (IP:Catalog No:112357, 4ug; Detection:Catalog No:112357 1:500) with HEK-293 cells lysate 1120ug.