

LTBP1 antibody

Catalog Number: 112359

Product name

LTBP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LTBP1 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human placenta tissue. Positive WB detected in mouse liver tissue, HeLa cells. Observed molecular weight by Western-blot: 150kd

Preparation

This antibody was obtained by immunization of LTBP1 recombinant protein (Accession Number: NM_000627). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:112359(LTBP1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112359(LTBP1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:112359(LTBP1 antibody) at dilution of 1:50

LTBP1 antibody



Catalog Number: 112359

(under 40x lens)