

LTBP2 antibody



Catalog Number: 112360

Product name

LTBP2 antibody

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

LTBP2 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, mouse testis tissue. Positive IHC detected in human placenta tissue. Observed molecular weight by Western-blot: 240 kDa




Preparation

This antibody was obtained by immunization of LTBP2 recombinant protein (Accession Number: NM_000428). Purification method: Antigen affinity purified.

human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112360(LTBP2 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human placenta using Catalog No:112360(LTBP2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:112360(LTBP2 antibody) at dilution of 1:50 (under 40x lens)