

# MAB21L2 antibody



Catalog Number: 112363

## Product name

Polyclonal

MAB21L2 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

MAB21L2 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 40 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of MAB21L2 recombinant protein (Accession Number: NM\_006439). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112363(MAB21L2 antibody) at dilution of 1:200