

MAdCAM1 antibody



Catalog Number: 112370

Product name

MAdCAM1 antibody

Specificity

Human; other species not tested.

Antibody description

MAdCAM1 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human small intestine tissue.

Preparation

This antibody was obtained by immunization of MAdCAM1 recombinant protein (Accession Number: BC142629). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:112370(MADCAM1 Antibody) at dilution of 1:200 (under 40x lens).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:112370(MADCAM1 Antibody) at dilution of 1:200 (under 40x lens).