

MAFB antibody



Catalog Number: 112373

Product name

Polyclonal

MAFB antibody

Ig Type

Specificity

Rabbit IgG

Human, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

MAFB Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 45 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of MAFB recombinant protein (Accession Number: NM_005461). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112373(MAFB Antibody) at dilution of 1:1000