

MAFB antibody



Catalog Number: 112373

Product name

MAFB antibody

Specificity

Human, Rat; other species not tested.

Antibody description

MAFB Rabbit Polyclonal antibody. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of MAFB recombinant protein (Accession Number: NM_005461). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112373(MAFB Antibody) at dilution of 1:1000