

MAGEA1 antibody



Catalog Number: 112376

Product name

MAGEA1 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEA1 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 43kd

Preparation

This antibody was obtained by immunization of MAGEA1 recombinant protein (Accession Number: NM_004988). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:112376(MAGEA1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:112376(MAGEA1 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:112376(MAGEA1 antibody) at dilution of 1:100 (under 40x lens)