

MAGEA1 antibody

Catalog Number: 112377

Product name

MAGEA1 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEA1 Rabbit Polyclonal antibody. Positive IF detected in A375 cells. Positive WB detected in A375 cells, MCF7 cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_004988). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:20-1:200

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:112377(MAGEA1 antibody) at dilution of 1:200



Immunofluorescent analysis of A375 cells using Catalog No:112377(MAGEA1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)