

MAGEB1 antibody

Catalog Number: 112384

Product name

MAGEB1 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEB1 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells, A375 cells. Positive WB detected in HepG2 cells, A375 cells. Observed molecular weight by Western-blot: 43-45 kDa

Preparation

This antibody was obtained by immunization of MAGEB1 recombinant protein (Accession Number: NM_002363). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112384(MAGEB1 antibody) at dilution of 1:2000



Immunofluorescent analysis of HepG2 cells using Catalog No:112384(MAGEB1 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)