

MAGEB4 antibody



Catalog Number: 112386

Product name

MAGEB4 antibody

Specificity

Human; other species not tested.

Antibody description

MAGEB4 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in K-562 cells. Observed molecular weight by Western-blot: 39 kDa

Preparation

This antibody was obtained by immunization of MAGEB4 recombinant protein (Accession Number: NM_002367). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:112386(MAGEB4 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using MAGEB4 antibody Catalog No:112386 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).