

# MAGED1 antibody



Catalog Number: 112389

## Product name

MAGED1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MAGED1 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, Apoptosised HeLa cells, HeLa cells, human brain tissue. Positive FC detected in HeLa cells. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 66 kDa, 92-95 kDa

## Preparation

This antibody was obtained by immunization of MAGED1 recombinant protein (Accession Number: NM\_001005333). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, FC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112389(MAGED1 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:112389(MAGED1 Antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemistry of paraffin-embedded human gliomas tissue slide using Catalog No:112389(MAGED1 Antibody) at dilution of

# MAGED1 antibody



Catalog Number: 112389

1:100 (under 40x lens)



1X10<sup>6</sup> HeLa cells were stained with 0.2ug MAGED1 antibody (Catalog No:112389, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.