

LUZP1 antibody

Catalog Number: 112395

Product name

LUZP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

LUZP1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 120 kDa

Preparation

This antibody was obtained by immunization of LUZP1 recombinant protein (Accession Number: BC051733). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112395(LUZP1 antibody) at dilution of 1:500



IP Result of anti-LUZP1 (IP:Catalog No:112395, 4ug; Detection:Catalog No:112395 1:500) with mouse brain tissue lysate 4000ug.