

# LXN antibody

Catalog Number: 112397

## Product name

LXN antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

LXN Rabbit Polyclonal antibody. Positive IP detected in mouse lung tissue. Positive WB detected in mouse lung tissue, human brain tissue, mouse ovary tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 26kd

## Preparation

This antibody was obtained by immunization of LXN recombinant protein (Accession Number: NM\_020169). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112397(LXN antibody) at dilution of 1:500



IP Result of anti-LXN (IP:Catalog No:112397, 4ug; Detection:Catalog No:112397 1:500) with mouse lung tissue lysate 5160ug.