

LY6K antibody

Catalog Number: 112403

Product name

LY6K antibody

Specificity

Human; other species not tested.

Antibody description

LY6K Rabbit Polyclonal antibody. Positive IP detected in PC-3 cells. Positive WB detected in human testis tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of LY6K recombinant protein (Accession Number: BC001291). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112403(LY6K antibody) at dilution of 1:500



IP Result of anti-LY6K (IP:Catalog No:112403, 4ug; Detection:Catalog No:112403 1:300) with PC-3 cells lysate 3200ug.