

LY96/MD2 antibody



Catalog Number: 112404

Product name

LY96/MD2 antibody

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

LY96/MD2 Rabbit Polyclonal antibody. Positive WB detected in U-937 cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 20 kDa and 30 kDa and 35 kDa

Preparation

This antibody was obtained by immunization of LY96/MD2 recombinant protein (Accession Number: NM_015364). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



U-937 cells were subjected to SDS PAGE followed by western blot with Catalog No:112404(LY96 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112404(LY96 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112404(LY96 antibody) at dilution of 1:50 (under 40x lens)