

# LYAR antibody



Catalog Number: 112405

## Product name

LYAR antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

LYAR Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HeLa cells, K-562 cells, mouse testis tissue. Observed molecular weight by Western-blot: 45 kDa

## Preparation

This antibody was obtained by immunization of LYAR recombinant protein (Accession Number: XM\_011513506). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112405(LYAR Antibody) at dilution of 1:1000



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:112405(LYAR Antibody) at dilution of 1:1000