

LYPLAL1 antibody



Catalog Number: 112411

Product name

LYPLAL1 antibody

WB: 1:200-1:2000

IHC: 1:20-1:200

Immunogen

[Human LYPLAL1 Recombinant protein \(GST tag\)](#)

Specificity

human, mouse; other species not tested.

Antibody description

LYPLAL1 Rabbit Polyclonal antibody. Positive IHC detected in human kidney. Positive WB detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 26kd

Preparation

This antibody was obtained by immunization of LYPLAL1 recombinant protein (Accession Number: XM_047445134). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112411(LYPLAL1 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112411(LYPLAL1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:112411(LYPLAL1

LYPLAL1 antibody



Catalog Number: 112411

antibody) at dilution of 1:50 (under 40x lens)