

LYSMD3 antibody



Catalog Number: 112413

Product name

LYSMD3 antibody

Immunogen

[Human LYSMD3 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

LYSMD3 Rabbit Polyclonal antibody. Positive WB detected in Neuro-2a cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of LYSMD3 recombinant protein (Accession Number: NM_198273). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:1000

Validations



Neuro-2a cells were subjected to SDS PAGE followed by western blot with Catalog No:112413(LYSMD3 Antibody) at dilution of 1:300