

LYVE1 antibody



Catalog Number: 112415

Product name

Polyclonal

LYVE1 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

LYVE1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 60 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of LYVE1 recombinant protein (Accession Number: NM_006691). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112415(LYVE1 antibody) at dilution of 1:500