

MAP4K3 antibody



Catalog Number: 112425

Product name

MAP4K3 antibody

Recommended Dilution:

WB: 1:1000-1:10000

Immunogen

IHC: 1:20-1:200

[Human MAP4K3 Recombinant protein \(His tag\)](#)


Validations

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAP4K3 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue, HeLa cells, human brain tissue. Positive IHC detected in human skeletal muscle tissue, human heart tissue. Observed molecular weight by Western-blot: 100kd



Preparation

This antibody was obtained by immunization of MAP4K3 recombinant protein (Accession Number: NM_001270425). Purification method: Antigen affinity purified.

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112425(MAP4K3 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Immunohistochemical of paraffin-embedded human skeletal muscle using Catalog No:112425(MAP4K3 antibody) at dilution of 1:100 (under 10x lens)



MAP4K3 antibody



Catalog Number: 112425

Immunohistochemical of paraffin-embedded
human skeletal muscle using Catalog

No:112425(MAP4K3 antibody) at dilution of 1:100
(under 40x lens)