

# MAPK7 antibody



Catalog Number: 112437

## Product name

MAPK7 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

MAPK7 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, human heart tissue, mouse brain tissue. Observed molecular weight by Western-blot: 59kd

## Preparation

This antibody was obtained by immunization of MAPK7 recombinant protein (Accession Number: XM\_011523957). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112437(MAPK7 antibody) at dilution of 1:300