

# MARCH8 antibody



Catalog Number: 112446

## Product name

MARCH8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MARCH8 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 30-33kd

## Preparation

This antibody was obtained by immunization of MARCH8 recombinant protein (Accession Number: XM\_047424767). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:112446(MARCH8 antibody) at dilution of



1:800

IP Result of anti-MARCH8 (IP:Catalog No:112446, 3ug; Detection:Catalog No:112446 1:400) with mouse lung tissue lysate 4000ug.