

Marcksl1 antibody



Catalog Number: 112450

Product name

Marcksl1 antibody

Specificity

human, mouse; other species not tested.

Antibody description

Marcksl1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of Recombinant Protein (Accession Number: NM_010807). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

158kd→

116kd→

97kd→

72kd→

56kd→

36kd→



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112450(Marcksl1 antibody) at dilution of 1:500