

# MAML3 antibody



Catalog Number: 112461

## Product name

MAML3 antibody

## Specificity

Human; other species not tested.

## Antibody description

MAML3 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 122 kDa

## Preparation

This antibody was obtained by immunization of MAML3 recombinant protein (Accession Number: NM\_018717). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:112461(MAML3 antibody) at dilution of 1:2000



Immunofluorescent analysis of HEK-293 cells using Catalog No:112461(MAML3 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)