

MAOA antibody



Catalog Number: 112469

Product name

WB: 1:200-1:2000

MAOA antibody

IP: 1:200-1:2000

Specificity

IHC: 1:20-1:200

Human; other species not tested.

Antibody description

MAOA Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue, human kidney tissue, human liver tissue, human pancreas tissue. Positive IP detected in human placenta tissue. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of MAOA recombinant protein (Accession Number: NM_000240). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

Recommended Dilution:


Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:112469(MAOA antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:112469(MAOA antibody) at dilution of 1:100 (under 25x lens)



MAOA antibody



Catalog Number: 112469

IP Result of anti-MAOA (IP:Catalog No:112469,

3ug; Detection:Catalog No:112469 1:600) with human placenta tissue lysate 2800ug.