

MARK3 antibody

Catalog Number: 112480

Product name

MARK3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MARK3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 80-90 kDa

Preparation

This antibody was obtained by immunization of MARK3 recombinant protein (Accession Number: NM_002376). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112480(CTAK1 antibody) at dilution of 1:200



IP Result of anti-CTAK1 (IP:Catalog No:112480, 3ug; Detection:Catalog No:112480 1:500) with mouse brain tissue lysate 4000ug.