

MARVELD2 antibody



Catalog Number: 112483

Product name

MARVELD2 antibody

Immunogen

[Human MARVELD2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MARVELD2 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells. Observed molecular weight by Western-blot: 64kd

Preparation

This antibody was obtained by immunization of MARVELD2 recombinant protein (Accession Number: XM_005248446). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:112483(MARVELD2 antibody) at dilution of 1:800