

MAST1 antibody



Catalog Number: 112490

Product name

MAST1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MAST1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 171 kDa

Preparation

This antibody was obtained by immunization of MAST1 recombinant protein (Accession Number: BC027985). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112490(MAST1 antibody) at dilution of 1:300



IP Result of anti-MAST1 (IP:Catalog No:112490, 4ug; Detection:Catalog No:112490 1:500) with mouse brain tissue lysate 4000ug.