

MAT2A antibody



Catalog Number: 112493

Product name

MAT2A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

MAT2A Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 40-45 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005911).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112493(MAT2A antibody) at dilution of 1:500



IP Result of anti-MAT2A (IP:Catalog No:112493, 4ug; Detection:Catalog No:112493 1:500) with HepG2 cells lysate 3200ug.