

# Matrilin 2 antibody



Catalog Number: 112498

---

## Product name

Matrilin 2 antibody

## Immunogen

[Human MATN2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Matrilin 2 Rabbit Polyclonal antibody. Positive WB detected in NIH/3T3 cells, HEK-293 cells, mouse liver tissue. Observed molecular weight by Western-blot: 120-150 kDa

## Preparation

This antibody was obtained by immunization of Matrilin 2 recombinant protein (Accession Number: NM\_030583). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:112498(MATN2 antibody) at dilution of 1:800