

MFAP5 antibody



Catalog Number: 112526

Product name

MFAP5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MFAP5 Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive WB detected in mouse lung tissue, Jurkat cells. Observed molecular weight by Western-blot: 20 kDa

Preparation

This antibody was obtained by immunization of MFAP5 recombinant protein (Accession Number: NM_003480). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:112526(MFAP5 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human placenta using Catalog No:112526(MFAP5 antibody) at dilution of 1:100 (under 40x lens)



Immunohistochemical of paraffin-embedded human placenta using Catalog No:112526(MFAP5 antibody) at dilution of 1:100 (under 10x lens)