

MFSD10 antibody



Catalog Number: 112531

Product name

MFSD10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MFSD10 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Observed molecular weight by Western-blot: 48-52 kDa

Preparation

This antibody was obtained by immunization of MFSD10 recombinant protein (Accession Number: NM_001146069). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:112531(MFSD10 Antibody) at dilution of 1:600