

MFSD7 antibody



Catalog Number: 112533

Product name

MFSD7 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

human, mouse, rat; other species not tested.

Applications

ELISA, WB

Antibody description

MFSD7 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells. Observed molecular weight by Western-blot: 46kd

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of MFSD7 recombinant protein (Accession Number: NM_032219). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:112533(MFSD7 antibody) at dilution of 1:400

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality