

# MBTD1 antibody



Catalog Number: 112540

## Product name

MBTD1 antibody

## Specificity

Human; other species not tested.

## Antibody description

MBTD1 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HL-60 cells, K-562 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 65-71 kDa

## Preparation

This antibody was obtained by immunization of MBTD1 recombinant protein (Accession Number: BC101736). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:112540(MBTD1 Antibody) at dilution of 1:600