

MC3R antibody

Catalog Number: 112541

Product name

MC3R antibody

Immunogen

[Human MC3R Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

MC3R Rabbit Polyclonal antibody. Positive IHC detected in human small intestine tissue.

Preparation

This antibody was obtained by immunization of MC3R recombinant protein (Accession Number: BC069599). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:112541(MC3R Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human small intestine tissue slide using Catalog No:112541(MC3R Antibody) at dilution of 1:50 (under 40x lens)