

MCHR1 antibody



Catalog Number: 112547

Product name

MCHR1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MCHR1 Rabbit Polyclonal antibody. Positive IHC detected in human heart tissue. Positive WB detected in SH-SY5Y cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of MCHR1 recombinant protein (Accession Number: BC021146). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:112547(MCHR1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112547(MCHR1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human heart tissue slide using Catalog No:112547(MCHR1 Antibody) at dilution of 1:50 (under 40x lens)