

# MCM4 antibody

Catalog Number: 112554

## Product name

MCM4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MCM4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, HL-60 cells, human liver tissue. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 97 kDa

## Preparation

This antibody was obtained by immunization of MCM4 recombinant protein (Accession Number: NM\_005914). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:112554(MCM4 antibody) at dilution of 1:600



IP Result of anti-MCM4 (IP:Catalog No:112554, 4ug; Detection:Catalog No:112554 1:600) with HEK-293 cells lysate 1200ug.