

Mecr antibody



Catalog Number: 112579

Product name

Mecr antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Mecr Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_025297). Purification method: Protein A purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:112579(Mecr antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:112579(Mecr antibody) at dilution of 1:50 (under 10x lens)